

Manufacturing Services

MQA provides custom manufacturing services for the biotechnology, pharmaceutical, medical device, and other life science industries. We are FDA registered and ISO 17025 accredited. MQA is available 7 days a week to provide you with a fast turnaround time.

Manufacturing Services

- Reagents
- Media
- Aseptic filling of injectables and IV bags for preclinical studies
- GLP & GMP cell bank preparation
- Aseptic assembly
- Steam sterilization
- Custom projects



GxP Compliance

- GxP compliance quality systems
- Validated clean room: ISO 5 BSC located in an ISO 7 cleanroom
- Validated equipment
- GMP trained personnel
- Qualification and media fill activities for aseptic filling follow FDA cGMP guidelines and are used to qualify aseptic process and personnel.
- Environmental controls, which include surface, viable, and non-viable particulate monitoring, the personnel who perform aseptic operations are glove sampled subsequent to completion of aseptic filling activities.
- QC testing



Microbiology & Quality Associates, Inc.
2341 Stanwell Drive, Concord, CA 94520
Office (925) 270-3800 Fax (925) 270-3806
Contact: sales@MicroQA.com
www.MicroQA.com